

---

# **QCPump**

*Release 0.1*

**Randle Taylor**

**Aug 01, 2022**



## CONTENTS:

<b>1</b>	<b>Release Notes</b>	<b>1</b>
1.1	v0.3.14 . . . . .	1
1.2	v0.3.13 . . . . .	1
1.3	v0.3.12 . . . . .	1
1.4	v0.3.11 . . . . .	2
1.5	v0.3.10 . . . . .	2
1.6	v0.3.9 . . . . .	2
1.7	v0.3.8 . . . . .	2
1.8	v0.3.6 . . . . .	2
1.9	v0.3.5 . . . . .	2
1.10	v0.3.4 . . . . .	3
1.11	v0.3.3 . . . . .	3
1.12	v0.3.2 . . . . .	3
1.13	v0.3.1 . . . . .	3
1.14	v0.3.0 . . . . .	3
<b>2</b>	<b>Installing QCPump</b>	<b>5</b>
2.1	Installing with the Windows Installer . . . . .	5
2.2	After Installing . . . . .	6
2.3	Starting QCPump Automatically . . . . .	6
2.4	Obtaining and running from source . . . . .	9
<b>3</b>	<b>QCPump Settings</b>	<b>11</b>
3.1	Available Settings . . . . .	12
<b>4</b>	<b>Operating QCPump</b>	<b>13</b>
4.1	Running your Pumps . . . . .	13
4.2	Stopping your Pumps . . . . .	13
4.3	Log Files, Config Files, and Settings Files . . . . .	13
<b>5</b>	<b>QCPump Built In Pump Types</b>	<b>15</b>
5.1	Configuring New Pumps . . . . .	15
5.2	Daily QA3 Pumps: Multiple Beams Per Test List . . . . .	36
5.3	QATrack+ Generic File Uploads . . . . .	45
<b>6</b>	<b>Development Notes</b>	<b>47</b>
6.1	Runing tests . . . . .	47
6.2	Release Checklist . . . . .	47
6.3	Building an exe on Windows . . . . .	48
6.4	Building the Installer . . . . .	48

<b>7</b>	<b>Pump Type Development</b>	<b>49</b>
7.1	Tutorial . . . . .	49
7.2	Developing your own Pump Types . . . . .	51
7.3	Adding Configuration Options To Your Pump . . . . .	52
7.4	Adding Validation To Your Pump . . . . .	52
7.5	Config Options . . . . .	55
7.6	Dependencies . . . . .	55
7.7	QATrack+ Mixins . . . . .	55
<b>8</b>	<b>QC Pump Overview</b>	<b>57</b>
8.1	QC Pump License . . . . .	58
<b>9</b>	<b>Indices and tables</b>	<b>59</b>

## RELEASE NOTES

### 1.1 v0.3.14

- Fix issue with making test level comments optional on MPC pumps

### 1.2 v0.3.13

- Made comment uploaded to QATrack+ optional (defaults to off because the comment prevents auto review)

### 1.3 v0.3.12

- Fixed bug related to logging of certificate paths
- Added a “Fast Search” option to MPC pumps. This option will restrict search for Results.csv files to subdirectories called MPCChecks. (defaults to on)
- Adjusted auth headers to make it possible for QCPump to talk to RadMachine
- For FFF beams, added Beam Shape Constancy results from the DQA3\_TREND table for DQA3 v1.06 SQL Server DBs. There are 9 new results sent to the server:

```
bsc{1,2,3}_{beam}  
bsc{1,2,3}_baseline_{beam}  
bsc{1,2,3}_diff_{beam}  
  
e.g.  
  
bsc2_6_fff  
bsc2_baseline_6_fff  
bsc_diff_6_fff
```

It is believed that bsc2 is the value shown in the DQA3 software so you should configure a test in QATrack+/RadMachine with slugs like bsc2\_6\_fff or bsc2\_10\_fff

- Partial work around for a suspected race condition occurring in the DQA3 database software. The race condition causes an incorrect data\_key to be written to the dqa3\_trend table, and hence data was being sent to the wrong unit. The workaround implemented here results in the dosimetry data being sent to the correct unit, but the temperature, pressure, comment, signature, and device serial number will still be incorrect if this occurs again. This is the result of a defect in the DQA3 software and should be exceedingly rare.

## 1.4 v0.3.11

- Fixed bug with base64 encoding of files for QATrackGenericBinaryFileUploader pump
- MPC pumps should now handle 2.5X and HDTSE beams
- MPC pumps should handle directory names like NDS-WKS-SN1234-2015-01-01-00-00-0001-10x-Beam
- QATrack+ validation requests receiving a 307 Temporary Redirect response will retry their request. This is an attempt (possibly in vain) to work around network monitoring software which may temporarily return 307s.

## 1.5 v0.3.10

- Added DQA SQL Server database pumps

## 1.6 v0.3.9

- Fix issue with missing Unit names for MPC
- Fix issue with “Â” appearing in some files
- Fixed broken link to QATrack docs

## 1.7 v0.3.8

- Allow disabling certificate patching by putting a file called nopatch.txt in C:\ProgramData\QATrack\Projectqcpump
- fix to retain unit choices when connection to QATrack+ fails for DQA3 pumps
- Attempt to fix SSL certificate errors in some networks.

## 1.8 v0.3.6

- Resolves an issue with unknown null urls being cached leading to pumps needing to be restarted in order to “re-discover” the URL to upload data to.

## 1.9 v0.3.5

- Ensure ssl certificate files are found / installed by pyinstaller
- Switch to using site\_config\_dir instead of user\_config\_dir for storing QCPump configs and logs so that multiple users can use same configuration. It is recommended that you “Install for All Users” when installing QCPump.
- Instructions for starting QCPump automatically have been added.

## 1.10 v0.3.4

- retracted

## 1.11 v0.3.3

- Fix stdio redirect
- Fix logging window history
- Add more json decoding errors

## 1.12 v0.3.2

- Handle non-missing test 400 Bad Requests when uploading to QATrack+

## 1.13 v0.3.1

- The File Mover pump types can now be set to *Copy* mode to have source files copied to the destination directory rather than moved
- The Python package *python-certifi-win32* has been added so that requests can use the Windows Certificate Store for SSL verification rather than using its bundled certificate chain. This should resolve some issues users were having with firewalls & network monitoring software.
- Added a missing permission for the DQA3 QCPump Firebird user

## 1.14 v0.3.0

### 1.14.1 General Features & Fixes

- A bug has been fixed where only the last subsection of a *Multiple* configuration section was being validated.
- Pumps which are marked as *Inactive* will not run validation code until they are re-activated. This eliminates un-necessary network calls and other validation work which, in addition to being more efficient overall, makes debugging QCPump itself simpler when multiple pumps are configured.
- New MPC pump type for uploading MPC results. See *Varian MPC Pumps*.
- New generic pump types for uploading Text & Binary Files to QATrack+ have been implemented. See *QATrack+ Generic File Uploads*.
- A *DISPLAY\_NAME* attribute has been added to Pump Types to aid with grouping together similar pump types when adding new pumps.
- Warning level debug messages were being logged as errors. This has been fixed.
- A new *PUMP\_ON\_STARTUP* (see *QCPump Settings*) setting has been added to allow pumping to begin immediately after QCPump is launched. This allows you to place QCPump in a startup folder and have it launched & start pumping when your computer is restarted.

### 1.14.2 DQA3 Pump Type Changes

- The *DATEADD* for calculating a *work\_completed* value in Firebird DQA3 queries has been eliminated in order to allow the query to work with Firebird versions < 2.1. *work\_completed* is now just calculated in Python code instead.
- The template for looking up Test Lists for beams now defaults to:

```
Daily QA3 Results: {{ beam_name }}
```

where *beam\_name* is the DQA3 test name (e.g. '6MeV', '6MV WDG', '6MV EDW 60 Weekly', '20 MeV DQA3 Daily'). This allows QCPump to handle a wider variety of beam types/configurations.

- More context variables are available when generating your test list name. In most cases you should only need to use *beam\_name*, however other variables are available should you need them. See the [DQA3 Test List Name docs](#).
- New *Multiple Beam Per Test List* DQA3 pumps have been added which will group results from multiple measurements together based on the results being recorded in a short window of time. There are two disadvantages to using the Multiple Beams Per Test List:
  1. If you have many beams configured this will result in long test lists which can impact performance when uploading data, or reviewing data in QATrack+.
  2. If you perform a measurement twice (e.g. take 2 6X measurements), only the 2nd result will be included.
- QATrack+ Unit names will now be displayed along with their Site in order to disambiguate units with the same name
- DQA3 machine names will now be shown with their Room name to disambiguate machines using the same tree names.



## INSTALLING QCPUMP

### 2.1 Installing with the Windows Installer

On Windows platforms please download the [QCPump Installer](#). Download and run the installer following the prompts. You may choose to install for the current user only, or all users on this system.




Fig. 1: QCPump: Choose whether to install for all users or just the current user.

Installing for All Users requires administrator privileges while installing for only the current user should be possible regardless of user type. Please note that either way you install, the pump configuration will be shared across all users regardless of whether you install for the current user or all users.

## 2.2 After Installing

Now that QCPump is installed, go to the *Configuring New Pumps* page to start configuring some *Pumps*. After that you may want to see the following section on starting QCPump automatically.

## 2.3 Starting QCPump Automatically

If you want want to have QCPump start up automatically when a user logs in you can do so by adding QCPump your startup folder.

First create a shortcut to QCPump on the desktop by right clicking on the desktop and selecting “New -> Shortcut”:




Fig. 2: QCPump: Create a desktop shortcut for QCPump

and select the QCPump application (usually located at *C:\Users\YOURUSERNAME\AppData\Local\Programs\QCPump\qcpump.exe* if you installed for a single user, or *C:\Program Files (x86)\QCPump\qcpump.exe* if you installed for all users).

and then finish the dialog.

Next we need to move the shortcut to the appropriate startup folder.

To startup for just the current user open the Run dialog (Win Key + R), type *shell:startup* and click OK.

Or to startup for all users open the Run dialog (Win Key + R), type *shell:common startup* and click OK.

Now drag and drop the desktop link into the startup folder:

Finally launch QCPump and ensure that “Run Pumps On Launch” is checked.




Fig. 3: Select the QCPump application to link to




Fig. 4: Open Startup folder




Fig. 5: Open Common Startup folder




Fig. 6: Startup Shortcut




Fig. 7: Tell QCPump to automatically start running pumps on launch

## 2.4 Obtaining and running from source

**Note:** You will need both Python (version 3.7-3.9) and Git installed on your computer to run QCPump from source.

In order to obtain the source code for QCPump install Git and then clone the QCPump repository:

```
git clone https://github.com/qatrackplus/qcpump.git
```

If you are on Windows, visit <https://www.lfd.uci.edu/~gohlke/pythonlibs/> and download DukPy for your particular version of Python. Now create a new venv to install the QCPump requirements:

```
cd qcpump

# create your venv
python -m venv env

# activate venv on Windows and install requirements
env\Scripts\Activate.ps1
pip install C:\path\to\dukpy-0.2.3-cp39-cp39-win_amd64.whl
pip install -r requirements\base.txt

# activate venv on *nix
source env/bin/activate
# replace 18.04 with your Ubuntu version
pip install -U -f https://extras.wxpython.org/wxPython4/extras/linux/gtk3/ubuntu-18.
↪04 wxPython
pip install -r requirements/base.txt
```

and then to run the program:


```
python launch_qcpump.py
```

Now you can proceed to the *Configuring New Pumps* page to start configuring some *Pumps*.




## QCPUMP SETTINGS

You can find the location of your QCPump settings file location file by launching QCPump and going to File->About:



You can edit the settings.json file with any text editor:



### 3.1 Available Settings

You may add one or more of the following settings to the settings.json file to override the default values:

**BROWSER\_USER\_AGENT (string)** The User-Agent to use when making outgoing web requests. (Default “Mozilla/5.0 (Windows NT 10.0; Win64; x64) AppleWebKit/537.36 (KHTML, like Gecko) Chrome/70.0.3538.102 Safari/537.36 Edge/18.19582”)

**DB\_CONNECT\_TIMEOUT (integer)** Timeout in seconds for database connections where available. (Default 30)

**DEBUG (true, false)** Currently only used to redirect std input / output. Must be set to *true* if you want to use an interactive debugger while developing QCPump. (Default: *false*)

**LOG\_LEVEL (debug, info, warning, error, critical)** Choose the QCPump application logging level (individual pumps are not affected by this). One of. (Default: *info*)

**LOG\_TO\_CONSOLE (true, false):** Should logs be written to console as well as log files?

**PUMP\_DIRECTORIES (list of file paths or null)** Set to list of other directories to include user defined pump types from. See *Pump Type Development*.

**PUMP\_ON\_STARTUP (true, false)** Should QCPump immediately start pumping when it is launched. This is useful for e.g. adding QCPump to a startup folder so it launches when a machine is rebooted and starts pumping immediately. (Default *false*)



## OPERATING QCPUMP

### 4.1 Running your Pumps

After you've configured and *saved at least one Pump* you can run all active pumps by clicking the *Run Pumps* toggle button at the top of the QCPump window.

#### 4.1.1 Running pumps automatically on launch

If you want your Pumps to start pumping automatically when they launch, set the *PUMP\_ON\_STARTUP* to *True*. This allows you to have QCPump launch automatically at Windows Startup and start pumping without user intervention.

### 4.2 Stopping your Pumps

At any time you may stop the current pumps by clicking the *Stop Pumps* toggle button at the top of the QCPump window.

### 4.3 Log Files, Config Files, and Settings Files

In order to see where QCPump is storing your log files, config files and QCPump setting files use the *About* menu option in the *File* menu:

Each *Pump* you configure will have its own configuration directory which contains a *config.json* JSON document that represents the current configuration of your *Pump*.

Each pump will have its own log file placed in the main QCPump log directory. These log files can provide important information when trying to debug any issues you are having with QCPump.

#### 4.3.1 QCPump Settings

QCPump has a limited set of settings you can configure by editing the *settings.json* JSON document found in the QCPump config directory. Settings available are:

**DEBUG** Useful for developers only. Set to "*DEBUG*": *true* to prevent redirecting stderr/stdout to file.

**LOG\_LEVEL** Controls the logging level for the QCPump application itself (not the pumps which have their own log level settings). Set to "*LOG\_LEVEL*": "*debug*" to get more detailed log info. Other options include *info*, *warning*, *error*, *critical*.

**DB\_CONNECTION\_TIMEOUT** Timeout for database connections where available. Default value is 3s.




Fig. 1: QCPump About Page

**PUMP\_DIRECTORIES** A list of directories to look for *custom Pump Types* you are using. This should be an empty list (`"PUMP_DIRECTORIES": []`) unless you are using custom *Pump Types* not included with QCPump.

A sample settings.json document might look like:

```
{
  "LOG_LEVEL": "info",
  "DEBUG": false,
  "PUMP_DIRECTORIES": ["C:/Users/yourusername/pumps/"],
  "DB_CONNECT_TIMEOUT": 3
}
```

## QCPUMP BUILT IN PUMP TYPES

See below for documentation relevant to specific pump types.

### 5.1 Configuring New Pumps

When you first launch QCPump you will be met with a more or less blank slate.




Fig. 1: QCPump with no pumps installed

Click the *Add Pump* button at the top right. Select the *Pump Type* you want to create, give the *Pump* a meaningful name and click *OK*.




Fig. 2: QCPump with no pumps installed

After adding you click *OK* you will see a tab for your new *Pump* displayed:

On the left of the tab are the configuration options for the pump and on the right is a logging window which is used to display messages from your pump.

The configuration panel is separated into different configuration sections. All pumps have a common set of configuration options included under the *Pump Configuration* heading which control things like whether this pump is active and how often the pump should run.

### 5.1.1 Common Pump Configuration Options

Every *Pump* you configure has the following options available in the *Pump Configuration* section:

**Type** Displays the *Pump Type* of the current *Pump* (not editable)

**Name** Displays the name of the current *Pump* (not editable)

**Interval (s)** How often the pump should run in seconds.

**Log Level** A dropdown allowing you set the verbosity of logging for this Pump. Set to *debug* to get the maximum level of verbosity. Logging messages will be shown both in the status panel as well as written to log files. The location of the log files can be seen by accessing the *About* menu option in the *File* menu.

**Active** Select whether this pump will be activated when you click the *Run Pumps* button.

*Pump Type* specific options are covered in their own pages:




Fig. 3: A new pump added

## FileMover Pumps

The two FileMover pumps are primarily meant to serve as simple examples of how to write your own pump types. There is both a *Simple File Mover* which simply moves all files from one directory to another, and a slightly more complex *File Mover* which allows you to filter based on file patterns, and iterate through subdirectories.

### Simple File Mover

The *Simple File Mover Pump* is a very simple pump that simply moves files from one or more source directories to one or more destination directories periodically. It will not recurse through subdirectories

**Warning:** Files in your Destination directory may be overwritten by new files if they have the same name!

### Configuration options

Define one or more file movers using these options:

**Source** The directory to move files from

**Destination:** The directory to move files to

**Mode:** Choose whether files should be moved, or copied




Fig. 4: QCPump Simple File Mover

## File Mover

The *File Mover Pump* is a pump that moves files from one or more source directories to one or more destination directories periodically. This pump differs from the *Simple File Mover* by adding options allowing you to perform *glob* file matching and filtering, as well as optionally recursing through subdirectories.

**Warning:** Files in your Destination directory may be overwritten by new files if they have the same name!

## Configuration options

Define one or more file movers using these options:

**Source** The directory to move files from

**Destination** The directory to move files to

**Mode:** Choose whether files should be moved, or copied

**Recursive** Set to *True* to have the file mover look through all subdirectories of the *Source* directory.

**Pattern** Set to a *Glob style pattern* to specify what files to include. Only files with a matching filename will be moved.

**Ignore Pattern** Set to a *Glob style pattern* to specify what files to exclude from the matched files. The *Ignore Pattern* is applied after the matching by *Pattern* is done and therefor acts as a filter on matched files. For example to move all *.dcm* files *except* for those that start with *foo-* set *Pattern* = *\*.dcm* and *Ignore Pattern* = *foo-*.




Fig. 5: QCPump File Mover

### Daily QA3 Pumps: One Beam Per Test List

This page refers to Pumps which upload results to a different test list for each beam. For example, if you have the following 10 beams configured on a single unit: 6X, 6FFF, 10X, 10FFF, 6X EDW60, 6E, 9E, 12E, 16E, 20E, then you will need 10 different test lists assigned to your unit in QATrack+ to record all of your DQA3 results. If you prefer to combine all of your results in a single test list, please see: *Daily QA3 Pumps: Multiple Beams Per Test List*.

---

**Note:** There are two disadvantages to using the Multiple Beams Per Test List:

1. If you have many beams configured this will result in very long test lists which can impact performance when uploading data, or reviewing data in QATrack+. 2. If you perform a measurement twice (e.g. take 2 6X measurements), only the 2nd result will be included.

---

QCPump currently has the ability to retrieve data from the following Daily QA3 data sources:

- DQA3 Firebird Database version 01.03
- DQA3 Firebird Database version 01.04
- DQA3 SQL Server Database version 01.06
- DQA3 data from Atlas 1.5

---

**Note:** The DQA3 pumps are tested on QATrack+ v3.1.X QCPump is not compatible with QATrack v0.3.X

---

#### Contents

- *Daily QA3 Pumps: One Beam Per Test List*
  - *Configuring QATrack+ for DQA3 Data*
  - *DQA3 Common Configuration Options*
  - *DQA3: Firebird Individual Beams Pump Type*
  - *DQA3: SQL Server Individual Beams DQA3 Pump Type*
  - *DQA3: Atlas Individual Beams DQA3 Pump Type*

### Configuring QATrack+ for DQA3 Data

In order to upload DQA3 data to QATrack+ you need to do a bit of setup work in QATrack+ first.

#### Create an API Token

In order to upload your data to QATrack+ via the API you will require an API token. See the [QATrack+ documentation](#) for how to create an API token. You may wish to create a dedicated user in QATrack+ just for use with QCPump. The user will only need a single permission in order to upload data: `qa | test list instance | Can add test list instance`.



## Configure Test Lists

**Note:** In order to simplify the creation of these test lists, there is a script included with QATrack+ v3.1.0+ to generate either a [Test Pack](#) or to create a TestList directly in your database. To run the script activate your virtualenv, changed to the QATrack+ root directory and then run

```
# create a 6X test list in the db (replace 6X with your beam type)
python manage.py runscript create_dqa3_testlist --script-args db 6X

# or create a 9E test pack (replace 9E with your beam type)
python manage.py runscript create_dqa3_testlist --script-args testpack 9E
```

QCPump requires QATrack+ to have a Test List configured for each beam type you want to upload results for. For example, if on your linacs you use 6X, 6FFF, 10X, 10FFF, 15X, 6E, 9E, 12E, 16E, 20E beams you will need 10 total test lists for DQA3 results. The Test List must have a specific set of attributes:

### Test List Name

The simplest method to have QCPump find the correct test list to upload data to is to give the test list a name which contains the name of the DQA3 Test (e.g. “6 MV”, or “6X Daily”, or “6 MV EDW60”) By default QCPump uses a test list name like:

Daily QA3 Results: {{ beam\_name }}

where QCPump will replace {{ beam\_name }} with the name of the DQA3 test beam (e.g. “Daily QA3 Results: 6X”) before searching QATrack+ for the test list to perform. You may also customize your test list name template with other variables which include:

**beam\_name** *dqa3\_template.tree\_name* from a Firebird database (e.g. “6 MV”, “6MeV”, “6MV WDG”, or *MachineTemplate.MachineTestName* from an Atlas database. e.g. (“6 MV DQA3 Daily”, “6MV EDW60 Weekly”, “20 MeV DQA3 Daily”)

**energy** The beam energy: 6, 9, 10, 15, 18 etc

**beam\_type** One of “Photon”, “FFF”, or “Electron”

**wedge\_type** Empty for non wedge beams, otherwise “EDW” or “Static”

**wedge\_angle** Empty for non wedge beams, otherwise 30, 45, 60 etc

**wedge\_orient** *dqa3\_template.wedgeorient* for FBD databases, or *MachineTemplate.WedgeOrient* for Atlas databases.

**device** Device serial number

**machine\_name** *dqa3\_machine.tree\_name* for FBD databases, or *Machine.MachineName* for Atlas databases.

**room\_name** *room.tree\_name* for FBD databases, or *Machine.RoomNumber* for Atlas databases.

In order to record your data in QATrack+ you will need to add tests to your Test list with one or more of the following macro names:

**data\_key: String** data\_key is a key from the DQA3 database used by QCPump and QATrack+ to ensure duplicate entries are not uploaded

**signature: String** signature is used to record the username of who completed the measurement

**temperature: Simple numerical** Temperature measured by the DQA3 device

**pressure: Simple numerical** The pressure measured by the DQA3 device

**dose:** **Simple Numerical** The dose measured by the DQA3 Device

**dose\_baseline:** **Simple Numerical** Baseline dose value used

**dose\_diff:** **Simple Numerical** Difference between measured dose and baseline

**axsym:** **Simple Numerical** Axial symmetry value

**axsym\_baseline:** **Simple Numerical** Axial symmetry baseline value

**axsym\_diff:** **Simple Numerical** Difference between measured axial symmetry and baseline

**trsym:** **Simple Numerical** Transverse symmetry value

**trsym\_baseline:** **Simple Numerical** Transverse symmetry baseline value

**trsym\_diff:** **Simple Numerical** Difference between measured transverse symmetry and baseline

**qaflat:** **Simple Numerical** Flatness value

**qaflat\_baseline:** **Simple Numerical** Flatness baseline value

**qaflat\_diff:** **Simple Numerical** Difference between measured flatness and baseline

**energy:** **Simple Numerical** Measured energy value

**energy\_baseline:** **Simple Numerical** Energy baseline value (always 0)

**energy\_diff:** **Simple Numerical** Difference between measured and baseline energy

**xsize:** **Simple Numerical** Measured width of profile in x direction

**xsize\_baseline:** **Simple Numerical** Baseline width of profile in x direction

**xsize\_diff:** **Simple Numerical** Difference between measured and baseline width of profile in x direction

**ysize:** **Simple Numerical** Measured width of profile in y direction

**ysize\_baseline:** **Simple Numerical** Baseline width of profile in y direction

**ysize\_diff:** **Simple Numerical** Difference between measured and baseline width of profile in y direction

**xshift:** **Simple Numerical** Measured shift of center of profile in x direction

**xshift\_baseline:** **Simple Numerical** Baseline shift of center of profile in x direction

**xshift\_diff:** **Simple Numerical** Difference between measured and baseline shift of center of profile in x direction

**yshift:** **Simple Numerical** Measured shift of center of profile in y direction

**yshift\_baseline:** **Simple Numerical** Baseline shift of center of profile in y direction

**yshift\_diff:** **Simple Numerical** Difference between measured and baseline shift of center of profile in y direction

Here is what a sample test list might look like:

### Assign Test Lists to Units

Once you have created these Test Lists in QATrack+ you need to [assign them to units](#) you want to record DQA3 data for.

**Name:**

**Slug:**   
A short unique name for use in the URL of this list

**Description:**   
A concise description of this test checklist. (You may use HTML markup)

**Description Preview:**

**Javascript:**   
Any extra javascript to run when loading perform page

**Warning message:**   
Message given when a test value is out of tolerance. Leave blank to disable warnings from being shown when tests are out of tolerance.

**Sublist Memberships**

This Test List is not a sublist of any other Test Lists.

TEST LIST MEMBERS		
TYPE	ID	MACRO NAME
Test	<input type="text" value="13888"/> Q x DQA3 Results: 6MV: 6MV: Signature	signature
Test	<input type="text" value="13889"/> Q x DQA3 Results: 6MV: 6MV: Temperature (°C)	temperature
Test	<input type="text" value="13890"/> Q x DQA3 Results: 6MV: 6MV: Pressure (kPa)	pressure
Test	<input type="text" value="13891"/> Q x DQA3 Results: 6MV: 6MV: Dose (%)	dose
Test	<input type="text" value="13892"/> Q x DQA3 Results: 6MV: 6MV: Dose Baseline (%)	dose_baseline
Test	<input type="text" value="13893"/> Q x DQA3 Results: 6MV: 6MV: Dose Diff (%)	dose_diff
Test	<input type="text" value="13894"/> Q x DQA3 Results: 6MV: 6MV: AxSym (%)	axsym
Test	<input type="text" value="13895"/> Q x DQA3 Results: 6MV: 6MV: AxSym Baseline (%)	axsym_baseline
Test	<input type="text" value="13896"/> Q x DQA3 Results: 6MV: 6MV: AxSym Diff (%)	axsym_diff
Test	<input type="text" value="13897"/> Q x DQA3 Results: 6MV: 6MV: TrSym (%)	trsym
Test	<input type="text" value="13898"/> Q x DQA3 Results: 6MV: 6MV: TrSym Baseline (%)	trsym_baseline
Test	<input type="text" value="13899"/> Q x DQA3 Results: 6MV: 6MV: TrSym Diff (%)	trsym_diff
Test	<input type="text" value="13900"/> Q x DQA3 Results: 6MV: 6MV: QAFIat (%)	qafiat
Test	<input type="text" value="13901"/> Q x DQA3 Results: 6MV: 6MV: QAFIat Baseline (%)	qafiat_baseline
Test	<input type="text" value="13902"/> Q x DQA3 Results: 6MV: 6MV: QAFIat Diff (%)	qafiat_diff

Fig. 6: A test list for recording 6MV results

## DQA3 Common Configuration Options

Most of the configuration options are the same for the two DQA3 *Pump Types*. Those settings are outlined here and the DQA3 database connection specific options are described below.

### QATrack+ API

**Api Url** Enter the root api url for the QATrack+ instance you want to upload data to. For Example <http://yourqatrackserver/api>

**Auth Token** Enter an authorization token for the QATrack+ instance you want to upload data to

**Throttle** Enter the minimum interval between data uploads (i.e. a value of 1 will allow 1 record per second to be uploaded)

**Verify SSL** Set to False if you want to bypass SSL certificate checks (e.g. if your QATrack+ instance is using a self signed certificate)

**Http Proxy** QCPump will try to autodetect your current proxy settings. However if you want to manually provide a proxy url you may do so. Proxy authentication url e.g. <http://10.10.1.10:3128> or <socks5://user:pass@host:port>

**Https Proxy** QCPump will try to autodetect your current proxy settings. However if you want to manually provide a proxy url you may do so. Proxy authentication url e.g. <https://10.10.1.10:3128> or <socks5://user:pass@host:port>

### Test List (depends on QATrack+ API)

**Name** Enter a template for searching QATrack+ for the name of the Test List you want to upload data to. The default is :

*Daily QA3 Results: {{ energy }}{{ beam\_type }}{{ wedge\_type }}{{ wedge\_angle }}*

In the template *{{ energy }}* will be replaced by the DQA3 beam energy (e.g. 6, 10, 15) and *{{ beam\_type }}* will be replaced by the DQA3 beam type (e.g. X, E, FFF). This template would result in QCPump trying to find a Test List called e.g. “Daily QA3 Results: 6X”.

### Unit (depends on QATrack+ API and DQA3Reader configs)

These config options are used to map DQA3 machine names to QATrack+ Unit names.

**Dqa3 Name** Select the DQA3 machine name to map

**Unit Name** Select the QATrack+ Unit name to map the DQA3 name to

### DQA3: Firebird Individual Beams Pump Type

Config options specific to Firebird DQA3 databases (01.03.00.00 & 01.04.00.00).

## DQA3Reader

**Host** Enter the host name of the Firebird database server you want to connect to

**Database** Enter the path to the database file you want to connect to on the server. For example  
C:\Users\YourUserName\databases\Snodata.fdb

**User** Enter the username you want to use to connect to the database with

**Password** Enter the password you want to use to connect to the database with

**Port** Enter the port number that the Firebird Database server is listening on

**Driver** Select the database driver you want to use. Use firebirdsql unless you have a good reason not to.

**History Days** Enter the number of prior days you want to look for data to import. If you are importing historical data you may want to temporarily set this to a large number of days (i.e. to get the last years worth of data set History days to 365) but normally a small number of days should be used to minimize the number of records fetched.

## Creating a Read-Only User for QCPump

While it is not required, you may wish to create a read only user for QCPump to connect to your database with. You may either use the Firebird tools *gsec* and *isql* to create the user or a third party tool like [FlameRobin](#) which is a great option for managing users and databases.

### Using gsec to create a new user

On the server where your Firebird database is located, open a CMD prompt and enter the following command to create a user with the username *qcpump* and password *qcpump*:

```
# for firebird 1.5
C:\Program Files (x86)\Firebird\Firebird_1_5\bin\gsec.exe" -user sysdba -password_
↵masterkey -database "localhost:C:\Program Files (x86)\Firebird\Firebird_1_5\
↵security.fdb

# for firebird 2.5
C:\Program Files (x86)\Firebird\Firebird_2_5\bin\gsec.exe" -user sysdba -password_
↵masterkey -database "localhost:C:\Program Files (x86)\Firebird\Firebird_1_5\
↵security2.fdb

GSEC> add qcpump -pw qcpump
GSEC> q
```

Next you can grant your user select rights using *isql*. Open *isql* specifying your username and password on the command line:

```
# for firebird 1.5
"C:\Program Files (x86)\Firebird\Firebird_1_5\bin\isql.exe" -user sysdba -password_
↵masterkey


# for firebird 2.5
"C:\Program Files (x86)\Firebird\Firebird_2_5\bin\isql.exe" -user sysdba -password_
↵masterkey
```

and connect to your database:

```
CONNECT "localhost:C:\Path\To\Your\Database\Snodata.fdb";
```

(note, you may need to replace `localhost`` with your actual server host name) then grant your user select rights on the tables required:

```
GRANT SELECT ON atlas_master to USER qcpump;
GRANT SELECT ON dqa3_machine to USER qcpump;
GRANT SELECT ON dqa3_trend to USER qcpump;
GRANT SELECT ON dqa3_data to USER qcpump;
GRANT SELECT ON device to USER qcpump;
GRANT SELECT ON dqa3_calibration to USER qcpump;
GRANT SELECT ON dqa3_template to USER qcpump;
GRANT SELECT ON dqa3_machine to USER qcpump;
GRANT SELECT ON room to USER qcpump;
quit;
```



```
Administrator: Command Prompt - "C:\Program Files (x86)\Firebird\Firebird_1_5\bin\isql.exe" -user sysdba -password masterkey
Microsoft Windows [Version 6.3.9600]
(c) 2013 Microsoft Corporation. All rights reserved.

C:\Users\Administrator>"C:\Program Files (x86)\Firebird\Firebird_1_5\bin\isql.exe" -user sysdba -password masterkey
Use CONNECT or CREATE DATABASE to specify a database
SQL> CONNECT "WIN-ALR5ISHQMJ6:C:\test.fdb";
Database: "WIN-ALR5ISHQMJ6:C:\test.fdb", User: sysdba
SQL> GRANT SELECT ON atlas_master to USER qcpump;
SQL>
```

Fig. 7: Grant qcpump user rights

You should now be able to use the username `qcpump` and password `qcpump` for the *User* and *Password* settings described above.

### DQA3: SQL Server Individual Beams DQA3 Pump Type

Config options specific to DQA3 SQL Server databases.

#### DQA3Reader

**Host** Enter the host name of the SQL Server database server you want to connect to

**Database** Enter the name of the database you want to connect to on the server. For example 'atlas'

**User** Enter the username you want to use to connect to the database with

**Password** Enter the password you want to use to connect to the database with

**Port** Enter the port number that the SQL Server database server is listening on

**Driver** Select the database driver you want to use. On Windows you will typically want to use the *ODBC Driver 17 for SQL Server* driver (ensure you have this driver installed on the computer running QCPump!). On Linux you will likely want to use one of the TDS drivers.

**History Days** Enter the number of prior days you want to look for data to import. If you are importing historical data you may want to temporarily set this to a large number of days (i.e. to get the last years worth of data set History days to 365) but normally a small number of days should be used to minimize the number of records fetched.

## DQA3: Atlas Individual Beams DQA3 Pump Type

Config options specific to Atlas DQA3 databases (SQLServer).

### DQA3Reader

**Host** Enter the host name of the SQL Server database server you want to connect to

**Database** Enter the name of the database you want to connect to on the server. For example 'atlas'

**User** Enter the username you want to use to connect to the database with

**Password** Enter the password you want to use to connect to the database with

**Port** Enter the port number that the SQL Server database server is listening on

**Driver** Select the database driver you want to use. On Windows you will typically want to use the *ODBC Driver 17 for SQL Server* driver (ensure you have this driver installed on the computer running QCPump!). On Linux you will likely want to use one of the TDS drivers.

**History Days** Enter the number of prior days you want to look for data to import. If you are importing historical data you may want to temporarily set this to a large number of days (i.e. to get the last years worth of data set History days to 365) but normally a small number of days should be used to minimize the number of records fetched.

## Varian MPC Pumps

### Contents

- *Varian MPC Pumps*
  - *How The MPC Pump Works*
  - *Configuring QATrack+ for MPC Data*
  - *MPC QCPump Configuration Options*

**Note:** The DQA3 pumps are tested on QATrack+ v3.1. QCPump is not compatible with QATrack v0.3.X

Varian's automated Machine Performance Check (MPC) stores its results in CSV files on disk which makes it easy to review those results using third party programs (thanks Varian!).

### How The MPC Pump Works

The MPC Pump works by periodically searching the *va\_transer\TDS* directory for any *Results.csv* files and then grouping the data from the csv files together and uploading that data to QATrack+. The data from the *Results.csv* files are combined based on the following three factors:

1. **The machine they were performed on.** QCPump uses the serial number from the MPCChecks directories to match with a QATrack+ Unit. For example, for an MPC directory like:

```
NDS-WKS-SN1234-2020-11-06-09-59-21-0009-GeometryCheckTemplate6xMVkV
```

QCPump considers the serial number to be 1234. In order to upload results to QATrack+, QCPump must be able to find a QATrack+ unit with its serial number set to 1234.

2. **The type of MPC check that was run.** QCPump will retrieve data from the following MPC directory types:

- a. **GeometryCheckTemplate6xMVkV** e.g. NDS-WKS-SN####-YYYY-MM-DD-HH-MM-SS-####-GeometryCheckTemplate6xMVkV\Results.csv
- b. **BeamCheckTemplate{energy}{beam\_type}** e.g. NDS-WKS-SN####-YYYY-MM-DD-HH-MM-SS-####-BeamCheckTemplate6x\Results.csv
- c. **GeometryCheckTemplate6xMVkVEnhancedCouch** e.g. NDS-WKS-SN####-YYYY-MM-DD-HH-MM-SS-####-GeometryCheckTemplate6xMVkVEnhancedCouch\Results.csv
- d. **EnhancedMLCCheckTemplate6x** e.g. NDS-WKS-SN####-YYYY-MM-DD-HH-MM-SS-####-EnhancedMLCCheckTemplate6x\Results.csv

results from a. and b. are grouped together into a single set of results to upload to QATrack+, while “enhanced” results from c. and d. are uploaded on their own.

3. **The date & time the results files were generated.** For the GeometryCheckTemplate & BeamCheckTemplate check types, QCPump will combine results from a period of N minutes into a single data set.

For example, given the following list of directories and Results.csv files with a time grouping window of 20 minutes:

```
NDS-WKS-SN5678-2020-06-25-07-10-30-0000-GeometryCheckTemplate6xMVkVEnhancedCouch\
↳Results.csv
NDS-WKS-SN5678-2020-06-25-07-20-30-0000-GeometryCheckTemplate6xMVkVEnhancedCouch\
↳Results.csv
NDS-WKS-SN6789-2020-09-30-09-15-17-0011-EnhancedMLCCheckTemplate6x\Results.csv
NDS-WKS-SN1234-2020-11-06-09-59-21-0009-GeometryCheckTemplate6xMVkV\Results.csv
NDS-WKS-SN1234-2020-11-06-09-59-21-0000-BeamCheckTemplate6xFFF\Results.csv
NDS-WKS-SN1234-2020-11-06-09-59-21-0000-BeamCheckTemplate10x\Results.csv
NDS-WKS-SN1234-2020-11-06-09-59-21-0000-BeamCheckTemplate10xFFF\Results.csv
NDS-WKS-SN1234-2020-11-06-10-01-21-0000-BeamCheckTemplate18x\Results.csv
NDS-WKS-SN1234-2020-11-06-10-01-21-0000-BeamCheckTemplate6e\Results.csv
NDS-WKS-SN1234-2020-11-06-10-01-21-0000-BeamCheckTemplate9e\Results.csv
NDS-WKS-SN1234-2020-11-06-10-02-21-0000-BeamCheckTemplate12e\Results.csv
NDS-WKS-SN1234-2020-11-06-10-02-21-0000-BeamCheckTemplate16e\Results.csv
NDS-WKS-SN1234-2020-11-06-10-02-21-0000-BeamCheckTemplate20e\Results.csv
NDS-WKS-SN1234-2020-11-06-10-03-21-0009-GeometryCheckTemplate6xMVkVEnhancedCouch\
↳Results.csv
NDS-WKS-SN1234-2020-11-06-11-59-21-0009-GeometryCheckTemplate6xMVkV\Results.csv
NDS-WKS-SN1234-2020-11-06-11-59-21-0000-BeamCheckTemplate6xFFF\Results.csv
NDS-WKS-SN1234-2020-11-06-11-59-21-0000-BeamCheckTemplate10x\Results.csv
NDS-WKS-SN1234-2020-11-06-11-59-21-0000-BeamCheckTemplate10xFFF\Results.csv
NDS-WKS-SN1234-2020-11-06-12-01-21-0000-BeamCheckTemplate18x\Results.csv
NDS-WKS-SN1234-2020-11-06-12-01-21-0000-BeamCheckTemplate6e\Results.csv
NDS-WKS-SN1234-2020-11-06-12-01-21-0000-BeamCheckTemplate9e\Results.csv
NDS-WKS-SN1234-2020-11-06-12-02-21-0000-BeamCheckTemplate12e\Results.csv
NDS-WKS-SN1234-2020-11-06-12-02-21-0000-BeamCheckTemplate16e\Results.csv
NDS-WKS-SN1234-2020-11-06-12-02-21-0000-BeamCheckTemplate20e\Results.csv
NDS-WKS-SN1234-2020-11-06-12-03-21-0009-GeometryCheckTemplate6xMVkVEnhancedCouch\
↳Results.csv
```

the results would be grouped as:

```
# Group 1) SN 5678, 2020-06-25 07:10, Enhanced Couch Checks
NDS-WKS-SN5678-2020-06-25-07-10-30-0000-GeometryCheckTemplate6xMVkVEnhancedCouch\
↳Results.csv
```

(continues on next page)



(continued from previous page)

```

# Group 2) SN 5678, 2020-06-25 07:20, Enhanced Couch Checks
# Even though these results were performed within 20min of Group 1),
# they are not grouped together because they are an "Enhanced" check type.
NDS-WKS-SN5678-2020-06-25-07-20-30-0000-GeometryCheckTemplate6xMVkVEnhancedCouch\
↳Results.csv

# Group 3) SN 6789, 2020-09-30 09:15, Enhanced MLC Checks
NDS-WKS-SN6789-2020-09-30-09-15-17-0011-EnhancedMLCCheckTemplate6x\Results.csv

# Group 4) SN 6789, 2020-11-06 09:59, Beam and Geometry Checks
# These are grouped together because they are not "enhanced" check types and
# all occurred within 20min of each other
NDS-WKS-SN1234-2020-11-06-09-59-21-0009-GeometryCheckTemplate6xMVkV\Results.csv
NDS-WKS-SN1234-2020-11-06-09-59-21-0000-BeamCheckTemplate6xFFF\Results.csv
NDS-WKS-SN1234-2020-11-06-09-59-21-0000-BeamCheckTemplate10x\Results.csv
NDS-WKS-SN1234-2020-11-06-09-59-21-0000-BeamCheckTemplate10xFFF\Results.csv
NDS-WKS-SN1234-2020-11-06-10-01-21-0000-BeamCheckTemplate18x\Results.csv
NDS-WKS-SN1234-2020-11-06-10-01-21-0000-BeamCheckTemplate6e\Results.csv
NDS-WKS-SN1234-2020-11-06-10-01-21-0000-BeamCheckTemplate9e\Results.csv
NDS-WKS-SN1234-2020-11-06-10-02-21-0000-BeamCheckTemplate12e\Results.csv
NDS-WKS-SN1234-2020-11-06-10-02-21-0000-BeamCheckTemplate16e\Results.csv
NDS-WKS-SN1234-2020-11-06-10-02-21-0000-BeamCheckTemplate20e\Results.csv

# Group 5) SN 1234, 2020-11-06 03:21, Enhanced Couch Checks
NDS-WKS-SN1234-2020-11-06-10-03-21-0009-GeometryCheckTemplate6xMVkVEnhancedCouch\
↳Results.csv

# Group 6) SN 1234, 2020-11-06 11:59, Beam and Geometry Checks
# These are grouped together because they are not "enhanced" check types and
# all occurred within 20min of each other. They are not grouped with Group 4)
# results because they occurred at least 20 min after the last result from Group 4
NDS-WKS-SN1234-2020-11-06-11-59-21-0009-GeometryCheckTemplate6xMVkV\Results.csv
NDS-WKS-SN1234-2020-11-06-11-59-21-0000-BeamCheckTemplate6xFFF\Results.csv
NDS-WKS-SN1234-2020-11-06-11-59-21-0000-BeamCheckTemplate10x\Results.csv
NDS-WKS-SN1234-2020-11-06-11-59-21-0000-BeamCheckTemplate10xFFF\Results.csv
NDS-WKS-SN1234-2020-11-06-12-01-21-0000-BeamCheckTemplate18x\Results.csv
NDS-WKS-SN1234-2020-11-06-12-01-21-0000-BeamCheckTemplate6e\Results.csv
NDS-WKS-SN1234-2020-11-06-12-01-21-0000-BeamCheckTemplate9e\Results.csv
NDS-WKS-SN1234-2020-11-06-12-02-21-0000-BeamCheckTemplate12e\Results.csv
NDS-WKS-SN1234-2020-11-06-12-02-21-0000-BeamCheckTemplate16e\Results.csv
NDS-WKS-SN1234-2020-11-06-12-02-21-0000-BeamCheckTemplate20e\Results.csv

# Group 7) SN 1234, 2020-11-06 12:03, Enhanced Couch Checks
NDS-WKS-SN1234-2020-11-06-12-03-21-0009-GeometryCheckTemplate6xMVkVEnhancedCouch\
↳Results.csv

```

## Configuring QATrack+ for MPC Data

In order to upload MPC data to QATrack+ you need to do a bit of setup work in QATrack+ first.

### Create an API Token

In order to upload your data to QATrack+ via the API you will require an API token. See the [QATrack+ documentation](#) for how to create an API token. You may wish to create a dedicated user in QATrack+ just for use with QCPump. The user will only need a single permission in order to upload data: *qa | test list instance | Can add test list instance*.

### Configure Your Unit's Serial Numbers

In order to determine which unit to upload MPC results to, QCPump queries the QATrack+ API to look for a unit with a serial number matching the MPC directory it finds results in. For example if QCPump finds a Results.csv file in a directory like “NDS-WKS-SN1234-2020-11-06-12-03-21-0009-GeometryCheckTemplate6xMVkVEnhancedCouch” then it will look for a unit configured in QATrack+ with a serial number of 1234. If QCPump can't find a unit with a matching serial number, the MPC results from that directory will be ignored.

### Configure Test Lists

You will need 1 or more test lists to record the MPC data uploaded by QCPump.

1. A test list named “MPC: Beam and Geometry Checks” to record data from *GeometryCheckTemplate6xMVkV & BeamCheckTemplate{energy}{beam\_type}* MPC results.
2. A test list named “MPC: Enhanced Couch Checks” to record data from *GeometryCheckTemplate6xMVkVEnhancedCouch* results files.
3. A test list named “MPC: Enhanced MLC Checks” to record data from *EnhancedMLCCheckTemplate6x* results files.

These test lists should have tests with macro names corresponding to the name of results in the Results.csv files. The Results.csv files have names like the following:

```
IsoCenterGroup/IsoCenterSize [mm]
IsoCenterGroup/IsoCenterMVOffset [mm]
IsoCenterGroup/IsoCenterKVOffset [mm]
BeamGroup/BeamOutputChange [%]
BeamGroup/BeamUniformityChange [%]
BeamGroup/BeamCenterShift [mm]
CollimationGroup/MLCGroup/MLCMaxOffsetA [mm]
CollimationGroup/MLCGroup/MLCMaxOffsetB [mm]
CollimationGroup/MLCGroup/MLCMeanOffsetA [mm]
GantryGroup/GantryAbsolute [°], -0.09, 0.3, Pass
GantryGroup/GantryRelative [°], 0.11, 0.3, Pass
```

and QCPump transforms these names into valid QATrack+ macro names using the following rules:

1. All slashes (/) and spaces are converted into underscores (\_)
2. Unit replacements are made as follows:
  - [mm] is replaced with *mm*,
  - [%] is replaced with *per*,

- [ $^{\circ}$ ] is replaced with *deg*
3. The beam energy/type is appended to the macro name
  4. The macro name is lowercased.

Examples of this substitution from a Results.csv file in a directory called *NDS-WKS-SN1234-2020-12-01-01-00-00-0009-GeometryCheckTemplate6xMVkV* are:

```
IsoCenterGroup/IsoCenterSize [mm] --> isocentergroup_isocentersize_mm_6x
BeamGroup/BeamOutputChange [%] --> beamgroup_beamoutputchange_per_6x
GantryGroup/GantryAbsolute [ $^{\circ}$ ] --> gantrygroup_gantryabsolute_deg_6x
```

**Warning:** Currently results for individual MLC leaves are not included. Any result which starts with any of these 4 strings:

```
CollimationGroup/MLCGroup/MLCLeavesA/MLCLeaf
CollimationGroup/MLCGroup/MLCLeavesB/MLCLeaf
CollimationGroup/MLCBacklashGroup/MLCBacklashLeavesA/MLCBacklashLeaf
CollimationGroup/MLCBacklashGroup/MLCBacklashLeavesB/MLCBacklashLeaf
```

will not be included in the api payload.

As a further example, the following Results.csv file found in a directory with the name *NDS-WKS-SN1234-2020-12-01-01-00-00-0009-GeometryCheckTemplate6xMVkV*

```
Name [Unit], Value, Threshold, Evaluation Result
IsoCenterGroup/IsoCenterSize [mm], 0.3, 0.5, Pass
IsoCenterGroup/IsoCenterMVOffset [mm], 0.26, 0.5, Pass
IsoCenterGroup/IsoCenterKVOffset [mm], 0.25, 0.5, Pass
BeamGroup/BeamOutputChange [%], -0.17, 2, Pass
BeamGroup/BeamUniformityChange [%], 0.31, 2, Pass
BeamGroup/BeamCenterShift [mm], 0.07, 0.5, Pass
CollimationGroup/MLCGroup/MLCMaxOffsetA [mm], 0.33, 1, Pass
CollimationGroup/MLCGroup/MLCMaxOffsetB [mm], 0.38, 1, Pass
CollimationGroup/MLCGroup/MLCMeanOffsetA [mm], 0.24, 1, Pass
...
GantryGroup/GantryAbsolute [ $^{\circ}$ ], -0.09, 0.3, Pass
```

would result in an API payload like this:

```
{
  "unit_test_collection": "https://qatrack.example.com/api/qa/unittestcollections/
↪1234/",
  "work_started": "2020-12-01 01:00",
  "work_completed": "2020-12-01 01:01",
  "user_key": "1234-2020-12-01-01-01",
  "day": 0,
  "tests": {
    isocentergroup_isocentersize_mm_6x: {"value": 0.3},
    isocentergroup_isocentermvoffset_mm_6x: {"value": 0.26},
    # ...
    beamgroup_beamoutputchange_per_6x: {"value": -0.17},
    # ... and so on
    gantrygroup_gantryabsolute_deg: {"value": -0.09}
    # ... and so on
  }
}
```

Therefore you will need to configure *Simple Numerical* tests for your test lists with these macro names (or a subset of them). The names of the tests can be anything you like, but naming your test the same as the names in the Results.csv file might be a good idea. So a Test List might look like:

Test List

ID: 806

Name:

Slug:   
A short unique name for use in the URL of this list

Description:   
A concise description of this test checklist. (You may use HTML markup)

Description Preview:

Javascript:   
Any extra javascript to run when loading perform page

Warning message:   
Message given when a test value is out of tolerance. Leave blank to disable warnings from being shown when tests are out of tolerance.

Sublist Memberships

This Test List is not a sublist of any other Test Lists.

TEST LIST MEMBERS

TYPE	ID		MACRO NAME	SHOW OUTLINE AROUND SUBLIST	REMOVE FROM TESTLIST
Test	<input type="text" value="12123"/>	IsoCenterGroup/IsoCenterSize [mm]: 6X	isocentergroup_isocentersize_mm_6x	-	<input type="checkbox"/>
Test	<input type="text" value="12124"/>	IsoCenterGroup/IsoCenterMVOffset [mm]: 6X	isocentergroup_isocentermvoffset_mm_6x	-	<input type="checkbox"/>
Test	<input type="text" value="12125"/>	IsoCenterGroup/IsoCenterKVOffset [mm]: 6X	isocentergroup_isocenterkvoffset_mm_6x	-	<input type="checkbox"/>
Test	<input type="text" value="12126"/>	BeamGroup/BeamOutputChange [%]: 6X	beamgroup_beamoutputchange_per_6x	-	<input type="checkbox"/>
Test	<input type="text" value="12127"/>	BeamGroup/BeamOutputChange [%]: 10X	beamgroup_beamoutputchange_per_10x	-	<input type="checkbox"/>
Test	<input type="text" value="12128"/>	BeamGroup/BeamOutputChange [%]: 18X	beamgroup_beamoutputchange_per_18x	-	<input type="checkbox"/>

Fig. 8: MPC Beam And Geometry Example Test List

**Note:** If QCPump detects that not all tests results for a given test list are included when it tries to upload results. It will automatically skip those results and attempt to upload the data again. This allows QCPump to upload partial result sets when e.g. you only run a single beam in MPC but your test list is configured to receive results from multiple beam types.

## Assign Test Lists to Units

Once you have created these Test Lists in QATrack+ you need to [assign them to units](#) you want to record DQA3 data for.

## MPC QCPump Configuration Options

### MPC

**TDS Directory** The “TDS directory” where MPC results are stored. Examples may be I:\TDS or \\YOURSERVER\VA\_Transer\TDS

**Fast Search** Restricts the search for Results.csv files to MPCChecks subdirectories. MPCChecks is the official Varian directory name and unless you have MPCResults in folders named something other than “MPCChecks”, it is recommended you leave this setting on for performance reasons.

**Days of history** The number of prior days you want to look for data to import. This should generally be 1 unless you are doing an initial import of historical results

**Results group time interval (min)** Enter the time interval (in minutes) for which results should be grouped together. That is to say, for Beam & Geometry checks how large of a time window should be used to consider MPC results part of the same session. This value should be a little bit longer than the typical time it takes you to run all your morning MPC checks.

**Wait for results (min)** Wait this many minutes for more results to be written to disk before uploading grouped results. In order to ensure all results from an MPC session, are written to disk, QCPump will wait this many minutes after the most recent Results.csv file it finds for a given machine before uploading results to QATrack+.

### QATrack+ API

**Api Url** Enter the root api url for the QATrack+ instance you want to upload data to. For Example <http://yourqatrackserver/api>

**Auth Token** Enter an authorization token for the QATrack+ instance you want to upload data to

**Throttle** Enter the minimum interval between data uploads (i.e. a value of 1 will allow 1 record per second to be uploaded)

**Verify SSL** Set to False if you want to bypass SSL certificate checks (e.g. if your QATrack+ instance is using a self signed certificate)

**Http Proxy** QCPump will try to autodetect your current proxy settings. However if you want to manually provide a proxy url you may do so. Proxy authentication url e.g. <http://10.10.1.10:3128> or socks5://user:pass@host:port

**Https Proxy** QCPump will try to autodetect your current proxy settings. However if you want to manually provide a proxy url you may do so. Proxy authentication url e.g. <https://10.10.1.10:3128> or socks5://user:pass@host:port

### Test List (depends on QATrack+ API)

**Name** Enter a template for searching QATrack+ for the name of the Test List you want to upload data to. The default is :

*MPC: {{ check\_type }}*

In the template *{{ check\_type }}* will be replaced by either:

- Beam and Geometry Checks
- Enhanced Couch Checks
- Enhanced MLC Checks

depending on the results being uploaded

### Configuration Validation

Any time you edit a configuration value your *Pump* will be re-validated. This may include actions such as ensuring a file system directory exists, or sending an HTTP request to check that a QATrack+ API endpoint is reachable.

---

**Note:** Currently the *Pump* tab is disabled during validation to prevent multiple validations from running simultaneously.

---

After validation, a message will be shown indicating whether each configuration setting is valid or not.

You can force a revalidation (say if some external factor was reconfigured / updated) by pressing the *Revalidate* button.

### Configuration Dependencies

Some *Pump Types* will have configuration sections which depend on other configuration sections being validated successfully before they can be completed. For example, selecting a Unit configured in QATrack+ will not be possible until QCPump has verified it can talk to your QATrack+ API.

When a dependency of a configuration section has not been met, the dependent section will be disabled and a message will be shown to inform you of which sections need to be completed before the dependent section can be verified.

### Multiple Configuration Subsections

Some configuration sections allow you to use multiple repeated subsections. An example of this would be the configuration used by the DQA3 Pump Types to map DQA3 Machines to QATrack+ Units:

To add new configuration subsections click the + button, and to remove a subsection, highlight one of its configuration options and click the -.

**QATrack+ API Configuration**

Validation: HTTPSPool(host='qatrack.example.com', port=443): Max retries exceeded with url: /api/auth/ (Caused by NewConnectionError('<urllib3.connection.HTTPSConnection object at 0x0000020D2E16A640>: Failed to establish a new connection: [Errno 11001] getaddrinfo failed'))

QATrack+ API	
Api Url	https://qatrack.example.com/api
Auth Token	97qf5be0362ef907dba79c24a31a6ndeobbnef54
Throttle	0.5
Verify Ssl	True
Http Proxy	
Https Proxy	

**Test List Configuration**

Please complete Config Section QATrack+ API to enable this section

Test List	
Name	Daily QA3 Results: {{ energy }}{{ beam_type }}
Data Key Test Name	DQA3 Results: Data Key

**Unit Configuration**

Please complete Config Sections QATrack+ API, DQA3Reader to enable this section

+ -

Unit #1	
Dqa3 Name	
QATrack+ Unit Name	

Fig. 9: QCPump config with unmet dependencies

+ -

Unit #1	
Dqa3 Name	
QATrack+ Unit Name	

Fig. 10: Multiple config subsections

## 5.1.2 Saving your Pumps

Before you can run your pumps, they all must be saved (pumps which haven't been saved will be marked with asterisks surrounding their pump tabs name). Once you have configured and validated your *Pump* click the *Save* button at the top of its config section. The location of the configuration files for QCPump can be found by accessing the *About* menu option in the *File* menu.

## 5.1.3 Resetting a Configuration

If you want to reset your Pump config, click the *Reset* button at the top of its config. This will reset the pumps config to the most recently saved version of the pumps configuration.

## 5.1.4 Deleting a Pump

If you want to delete a Pump config, click the *Delete* button at the top of its config. This will permanently delete this pump and remove its configuration file from disk.

## 5.2 Daily QA3 Pumps: Multiple Beams Per Test List

This page refers to Pumps which group together results from a single unit recorded during a short time period to a single test list. For example, if you have the following 10 beams configured on a single unit: 6X, 6FFF, 10X, 10FFF, 6X EDW60, 6E, 9E, 12E, 16E, 20E, and perform measurements for all these beams in a 10 minute period, then QCPump will combine all of these results and upload to QATrack+ to perform a single test list. If you prefer to use a single test list per beam type to record your results, please see: *Daily QA3 Pumps: One Beam Per Test List*.

---

**Note:** There are two disadvantages to using the Multiple Beams Per Test List:

1. If you have many beams configured this will result in very long test lists which can impact performance when uploading data, or reviewing data in QATrack+.
  2. If you perform a measurement twice (e.g. take 2 6X measurements), only the 2nd result will be included.
- 

QCPump currently has the ability to retrieve data from the following Daily QA3 data sources:

- DQA3 Firebird Database version 01.03
  - DQA3 Firebird Database version 01.04
  - DQA3 SQL Server Database version 01.06
  - DQA3 data from Atlas 1.5
- 

**Note:** The DQA3 pumps are tested on QATrack+ v3.1.X QCPump is not compatible with QATrack v0.3.X

---

### Contents

- *Daily QA3 Pumps: Multiple Beams Per Test List*
  - *Configuring QATrack+ for DQA3 Data*
  - *DQA3 Common Configuration Options*



- *DQA3: Firebird Multiple Beams Per Test List Pump Type*
- *DQA3: SQL Server Multiple Beams Per Test List Pump Type*
- *DQA3: Atlas Multiple Beams Per Test List Pump Type*

## 5.2.1 Configuring QATrack+ for DQA3 Data

In order to upload DQA3 data to QATrack+ you need to do a bit of setup work in QATrack+ first.

### Create an API Token

In order to upload your data to QATrack+ via the API you will require an API token. See the [QATrack+ documentation](#) for how to create an API token. You may wish to create a dedicated user in QATrack+ just for use with QCPump. The user will only need a single permission in order to upload data: *qa | test list instance | Can add test list instance*.

### Configure Test Lists

**Note:** In order to simplify the creation of these test lists, there is a script included with QATrack+ v3.1.0+ to generate either a [Test Pack](#) or to create a TestList directly in your database. To run the script activate your virtualenv, changed to the QATrack+ root directory and then run

```
# create a test list in the db for photons
python manage.py runscript create_grouped_dqa3_testlist --script-args db "Daily QA3_
↳Results: Photons" 6X 6FFF 10X 10FFF "6X EDW60"

# or create a test pack
python manage.py runscript create_grouped_dqa3_testlist --script-args testpack "Daily_
↳QA3 Results: Electrons" 6E 9E 12E 16E 20E
```

QCPump requires QATrack+ to have a single Test List per unit configured to record all your Daily QA3 Results. The Test List must have a specific set of attributes:

**Test List Name** Enter the name of the test list you configured in QATrack+ to record your DQA3 results e.g.:

Daily QA3 Results: Photons

When QCPump gathers data to post to QATrack+, it will convert the DQA3 test name, into a QATrack+ macro name for each test result. For example if the DQA3 Test/Beam is named “6 MV EDW60” then QCPump would generate test results for macros named *dose\_6\_mv\_edw60*, *dose\_baseline\_6\_mv\_edw60*, *dose\_diff\_6\_mv\_edw60* and so on. You will need to create tests with at least one of the following macro names (where {beam\_name} is suitably replaced according to which beam it is for):

**data\_key\_{beam\_name}: String** data\_key is a key from the DQA3 database used by QCPump and QATrack+ to ensure duplicate entries are not uploaded

**signature\_{beam\_name}: String** signature is used to record the username of who completed the measurement

**temperature\_{beam\_name}: Simple numerical** Temperature measured by the DQA3 device

**pressure\_{beam\_name}: Simple numerical** The pressure measured by the DQA3 device

**dose\_{beam\_name}: Simple Numerical** The dose measured by the DQA3 Device

**dose\_baseline\_{beam\_name}: Simple Numerical** Baseline dose value used

**dose\_diff\_{beam\_name}: Simple Numerical** Difference between measured dose and baseline

**axsym\_{beam\_name}: Simple Numerical** Axial symmetry value

**axsym\_baseline\_{beam\_name}: Simple Numerical** Axial symmetry baseline value

**axsym\_diff\_{beam\_name}: Simple Numerical** Difference between measured axial symmetry and baseline

**trsym\_{beam\_name}: Simple Numerical** Transverse symmetry value

**trsym\_baseline\_{beam\_name}: Simple Numerical** Transverse symmetry baseline value

**trsym\_diff\_{beam\_name}: Simple Numerical** Difference between measured transverse symmetry and baseline

**qaflat\_{beam\_name}: Simple Numerical** Flatness value

**qaflat\_baseline\_{beam\_name}: Simple Numerical** Flatness baseline value

**qaflat\_diff\_{beam\_name}: Simple Numerical** Difference between measured flatness and baseline

**energy\_{beam\_name}: Simple Numerical** Measured energy value

**energy\_baseline\_{beam\_name}: Simple Numerical** Energy baseline value (always 0)

**energy\_diff\_{beam\_name}: Simple Numerical** Difference between measured and baseline energy

**xsize\_{beam\_name}: Simple Numerical** Measured width of profile in x direction

**xsize\_baseline\_{beam\_name}: Simple Numerical** Baseline width of profile in x direction

**xsize\_diff\_{beam\_name}: Simple Numerical** Difference between measured and baseline width of profile in x direction

**ysize\_{beam\_name}: Simple Numerical** Measured width of profile in y direction

**ysize\_baseline\_{beam\_name}: Simple Numerical** Baseline width of profile in y direction

**ysize\_diff\_{beam\_name}: Simple Numerical** Difference between measured and baseline width of profile in y direction

**xshift\_{beam\_name}: Simple Numerical** Measured shift of center of profile in x direction

**xshift\_baseline\_{beam\_name}: Simple Numerical** Baseline shift of center of profile in x direction

**xshift\_diff\_{beam\_name}: Simple Numerical** Difference between measured and baseline shift of center of profile in x direction


**yshift\_{beam\_name}: Simple Numerical** Measured shift of center of profile in y direction

**yshift\_baseline\_{beam\_name}: Simple Numerical** Baseline shift of center of profile in y direction

**yshift\_diff\_{beam\_name}: Simple Numerical** Difference between measured and baseline shift of center of profile in y direction

Here is an example of what a test list configured with a sublist per beam might look like:

and the sublist for recording the 6MV results:

**Name:**  

---

**Slug:**   
A short unique name for use in the URL of this list

---

**Description:**   
A concise description of this test checklist. (You may use HTML markup)

---

**Description Preview:**

---

**Javascript:**   
Any extra javascript to run when loading perform page

---

**Warning message:**   
Message given when a test value is out of tolerance. Leave blank to disable warnings from

**Sublist Memberships**

This Test List is not a sublist of any other Test Lists.













TEST LIST MEMBERS		
TYPE	ID	
Sublist	<input type="text" value="838"/>	Q x (838) Daily QA3 Results: 6MV Daily  
Sublist	<input type="text" value="839"/>	Q x (839) Daily QA3 Results: 10MV Daily  
Sublist	<input type="text" value="840"/>	Q x (840) Daily QA3 Results: 10FFF Daily  
Sublist	<input type="text" value="841"/>	Q x (841) Daily QA3 Results: 6MV - 15EDW  
Sublist	<input type="text" value="842"/>	Q x (842) Daily QA3 Results: 6MV - IMRT daily  
Sublist	<input type="text" value="843"/>	Q x (843) Daily QA3 Results: 6MV - 60EDW  
Test	<input type="text"/>	

Fig. 11: Parent test list for recording DQA3 results

**Name:**

**Slug:**   
A short unique name for use in the URL of this list

**Description:**   
A concise description of this test checklist. (You may use HTML markup)

**Description Preview:**

**Javascript:**   
Any extra javascript to run when loading perform page

**Warning message:**   
Message given when a test value is out of tolerance. Leave blank to disable warnings from being shown when tests are out of tolerance.

**Sublist Memberships**

**This Test List is a Sublist of the following Test Lists:**

- Daily QA3 Results: Atlas

TEST LIST MEMBERS		
TYPE	ID	MACRO NAME
Test	<input type="text" value="13702"/> Q x DQA3 Results: 6MV Daily: Signature ✓ +	signature_6mv_daily
Test	<input type="text" value="13703"/> Q x DQA3 Results: 6MV Daily: Temperature (°C) ✓ +	temperature_6mv_daily
Test	<input type="text" value="13704"/> Q x DQA3 Results: 6MV Daily: Pressure (kPa) ✓ +	pressure_6mv_daily
Test	<input type="text" value="13705"/> Q x DQA3 Results: 6MV Daily: Dose (%) ✓ +	dose_6mv_daily
Test	<input type="text" value="13706"/> Q x DQA3 Results: 6MV Daily: Dose Baseline (%) ✓ +	dose_baseline_6mv_daily
Test	<input type="text" value="13707"/> Q x DQA3 Results: 6MV Daily: Dose Diff (%) ✓ +	dose_diff_6mv_daily
Test	<input type="text" value="13708"/> Q x DQA3 Results: 6MV Daily: AxSym (%) ✓ +	axsym_6mv_daily
Test	<input type="text" value="13709"/> Q x DQA3 Results: 6MV Daily: AxSym Baseline (%) ✓ +	axsym_baseline_6mv_daily
Test	<input type="text" value="13710"/> Q x DQA3 Results: 6MV Daily: AxSym Diff (%) ✓ +	axsym_diff_6mv_daily
Test	<input type="text" value="13711"/> Q x DQA3 Results: 6MV Daily: TrSym (%) ✓ +	trsym_6mv_daily
Test	<input type="text" value="13712"/> Q x DQA3 Results: 6MV Daily: TrSym Baseline (%) ✓ +	trsym_baseline_6mv_daily
Test	<input type="text" value="13713"/> Q x DQA3 Results: 6MV Daily: TrSym Diff (%) ✓ +	trsym_diff_6mv_daily
Test	<input type="text" value="13714"/> Q x DQA3 Results: 6MV Daily: OAFlat (%) ✓ +	oaflat_6mv_daily

Fig. 12: Child test list for recording 6MV DQA3 results

## Assign Test Lists to Units

Once you have created these Test Lists in QATrack+ you need to **assign them to units** you want to record DQA3 data for.

### 5.2.2 DQA3 Common Configuration Options

Most of the configuration options are the same for the two DQA3 *Pump Types*. Those settings are outlined here and the DQA3 database connection specific options are described below.

#### QATrack+ API

**Api Url** Enter the root api url for the QATrack+ instance you want to upload data to. For Example `http://yourqatrackserver/api`

**Auth Token** Enter an authorization token for the QATrack+ instance you want to upload data to

**Throttle** Enter the minimum interval between data uploads (i.e. a value of 1 will allow 1 record per second to be uplodged)

**Verify SSL** Set to False if you want to bypass SSL certificate checks (e.g. if your QATrack+ instance is using a self signed certificate)

**Http Proxy** QCPump will try to autodetect your current proxy settings. However if you want to manually provide a proxy url you may do so. Proxy authentication url e.g. `http://10.10.1.10:3128` or `socks5://user:pass@host:port`

**Https Proxy** QCPump will try to autodetect your current proxy settings. However if you want to manually provide a proxy url you may do so. Proxy authentication url e.g. `https://10.10.1.10:3128` or `socks5://user:pass@host:port`

#### Test List (depends on QATrack+ API)

**Name** Enter a template for searching QATrack+ for the name of the Test List you want to upload data to. The default is :

*Daily QA3 Results: {{ energy }}{{ beam\_type }}*

In the template `{{ energy }}` will be replaced by the DQA3 beam energy (e.g. 6, 10, 15) and `{{ beam_type }}` will be replaced by the DQA3 beam type (e.g. X, E, FFF). This template would result in QCPump trying to find a Test List called e.g. "Daily QA3 Results: 6X".

#### Unit (depends on QATrack+ API and DQA3Reader configs)

These config options are used to map DQA3 machine names to QATrack+ Unit names.

**Dqa3 Name** Select the DQA3 machine name to map

**Unit Name** Select the QATrack+ Unit name to map the DQA3 name to

### 5.2.3 DQA3: Firebird Multiple Beams Per Test List Pump Type

Config options specific to Firebird DQA3 databases (01.03.00.00 & 01.04.00.00).

#### DQA3Reader

**Host** Enter the host name of the Firebird database server you want to connect to

**Database** Enter the path to the database file you want to connect to on the server. For example  
C:\Users\YourUserName\databases\Sncdata.fdb

**User** Enter the username you want to use to connect to the database with

**Password** Enter the password you want to use to connect to the database with

**Port** Enter the port number that the Firebird Database server is listening on

**Driver** Select the database driver you want to use. Use firebirdsql unless you have a good reason not to.

**History Days** Enter the number of prior days you want to look for data to import. If you are importing historical data you may want to temporarily set this to a large number of days (i.e. to get the last years worth of data set History days to 365) but normally a small number of days should be used to minimize the number of records fetched.

**Results group time interval (min)** Enter the time interval (in minutes) for which results should be grouped together. That is to say, for Beam & Geometry checks how large of a time window should be used to consider MPC results part of the same session. This value should be a little bit longer than the typical time it takes you to run all

**Wait for results (min)** Wait this many minutes for more results to be written to disk before uploading grouped results. In order to ensure all results from an MPC session, are written to disk, QCPump will wait this many minutes after the most recent Results.csv file it finds for a given machine before uploading results to QATrack+.

**Beam Types** Select which beam types (Photon, Electron, All) you want to include for this pump. It is a good idea to create separate pumps for electrons and photons and corresponding test lists in QATrack+ for recording photon & electron results separately. The Photon option will include FFF & wedged beams.

#### Creating a Read-Only User for QCPump

While it is not required, you may wish to create a read only user for QCPump to connect to your database with. You may either use the Firebird tools *gsec* and *isql* to create the user or a third party tool like [FlameRobin](#) which is a great option for managing users and databases.

##### Using gsec to create a new user

On the server where your Firebird database is located, open a CMD prompt and enter the following command to create a user with the username *qcpump* and password *qcpump*:

```
# for firebird 1.5
C:\Program Files (x86)\Firebird\Firebird_1_5\bin\gsec.exe" -user sysdba -password_
↪masterkey -database "localhost:C:\Program Files (x86)\Firebird\Firebird_1_5\
↪security.fdb

# for firebird 2.5
C:\Program Files (x86)\Firebird\Firebird_2_5\bin\gsec.exe" -user sysdba -password_
↪masterkey -database "localhost:C:\Program Files (x86)\Firebird\Firebird_1_5\
↪security2.fdb

GSEC> add qcpump -pw qcpump
GSEC> q
```

Next you can grant your user select rights using isql. Open isql specifying your username and password on the command line:

```
# for firebird 1.5
"C:\Program Files (x86)\Firebird\Firebird_1_5\bin\isql.exe" -user sysdba -password_
↵masterkey

# for firebird 2.5
"C:\Program Files (x86)\Firebird\Firebird_2_5\bin\isql.exe" -user sysdba -password_
↵masterkey
```

and connect to your database:

```
CONNECT "localhost:C:\Path\To\Your\Database\Srcdata.fdb";
```

(note, you may need to replace `localhost` with your actual server host name) then grant your user select rights on the tables required:

```
GRANT SELECT ON atlas_master to USER qcpump;
GRANT SELECT ON dqa3_machine to USER qcpump;
GRANT SELECT ON dqa3_trend to USER qcpump;
GRANT SELECT ON dqa3_data to USER qcpump;
GRANT SELECT ON device to USER qcpump;
GRANT SELECT ON dqa3_calibration to USER qcpump;
GRANT SELECT ON dqa3_template to USER qcpump;
GRANT SELECT ON dqa3_machine to USER qcpump;
GRANT SELECT ON room to USER qcpump;
quit;
```

```
Administrator: Command Prompt - "C:\Program Files (x86)\Firebird\Firebird_1_5\bin\isql.exe" -user sysdba -password masterkey
Microsoft Windows [Version 6.3.9600]
(c) 2013 Microsoft Corporation. All rights reserved.

C:\Users\Administrator>"C:\Program Files (x86)\Firebird\Firebird_1_5\bin\isql.exe" -user sysdba -password masterkey
Use CONNECT or CREATE DATABASE to specify a database
SQL> CONNECT "WIN-ALR51SHQM16\C:\test.fdb";
Database: "WIN-ALR51SHQM16\C:\test.fdb", User: sysdba
SQL> GRANT SELECT ON atlas_master to USER qcpump;
SQL>
```

Fig. 13: Grant qcpump user rights

You should now be able to use the username `qcpump` and password `qcpump` for the *User* and *Password* settings described above.

## 5.2.4 DQA3: SQL Server Multiple Beams Per Test List Pump Type

Config options specific to SQL Server DQA3 databases.

### DQA3Reader

**Host** Enter the host name of the SQL Server database server you want to connect to

**Database** Enter the name of the database you want to connect to on the server. For example ‘atlas’

**User** Enter the username you want to use to connect to the database with

**Password** Enter the password you want to use to connect to the database with

**Port** Enter the port number that the SQL Server database server is listening on

**Driver** Select the database driver you want to use. On Windows you will typically want to use the *ODBC Driver 17 for SQL Server* driver (ensure you have this driver installed on the computer running QCPump!). On Linux you will likely want to use one of the TDS drivers.

**History Days** Enter the number of prior days you want to look for data to import. If you are importing historical data you may want to temporarily set this to a large number of days (i.e. to get the last years worth of data set History days to 365) but normally a small number of days should be used to minimize the number of records fetched.

**Results group time interval (min)** Enter the time interval (in minutes) for which results should be grouped together. That is to say, for Beam & Geometry checks how large of a time window should be used to consider MPC results part of the same session. This value should be a little bit longer than the typical time it takes you to run all

**Wait for results (min)** Wait this many minutes for more results to be written to disk before uploading grouped results. In order to ensure all results from an MPC session, are written to disk, QCPump will wait this many minutes after the most recent Results.csv file it finds for a given machine before uploading results to QATrack+.

**Beam Types** Select which beam types (Photon, Electron, All) you want to include for this pump. It is a good idea to create separate pumps for electrons and photons and corresponding test lists in QATrack+ for recording photon & electron results seperately. The Photon option will include FFF & wedged beams.

## 5.2.5 DQA3: Atlas Multiple Beams Per Test List Pump Type

Config options specific to Atlas DQA3 databases (SQLServer).

### DQA3Reader

**Host** Enter the host name of the SQL Server database server you want to connect to

**Database** Enter the name of the database you want to connect to on the server. For example 'atlas'

**User** Enter the username you want to use to connect to the database with

**Password** Enter the password you want to use to connect to the database with

**Port** Enter the port number that the SQL Server database server is listening on

**Driver** Select the database driver you want to use. On Windows you will typically want to use the *ODBC Driver 17 for SQL Server* driver (ensure you have this driver installed on the computer running QCPump!). On Linux you will likely want to use one of the TDS drivers.

**History Days** Enter the number of prior days you want to look for data to import. If you are importing historical data you may want to temporarily set this to a large number of days (i.e. to get the last years worth of data set History days to 365) but normally a small number of days should be used to minimize the number of records fetched.

**Results group time interval (min)** Enter the time interval (in minutes) for which results should be grouped together. That is to say, for Beam & Geometry checks how large of a time window should be used to consider MPC results part of the same session. This value should be a little bit longer than the typical time it takes you to run all

**Wait for results (min)** Wait this many minutes for more results to be written to disk before uploading grouped results. In order to ensure all results from an MPC session, are written to disk, QCPump will wait this many minutes after the most recent Results.csv file it finds for a given machine before uploading results to QATrack+.

**Beam Types** Select which beam types (Photon, Electron, All) you want to include for this pump. It is a good idea to create separate pumps for electrons and photons and corresponding test lists in QATrack+ for recording photon & electron results seperately. The Photon option will include FFF & wedged beams.



## 5.3 QATrack+ Generic File Uploads

QCPump currently has two pumps for uploading text or binary files to QATrack+. Both pumps operate by watching a directory for files, uploading them to a QATrack+ test list, and optionally, moving the file to a new directory after the file is processed.

### 5.3.1 Configuration options

---

**Todo:** Add common config options

---

<p><b>Warning:</b> Files in your Destination directory may be overwritten by new files if they have the same name!</p>
--

#### QATrack+ File Upload: Generic Text File

For uploading text files (e.g. csv, Profiler exports, etc)

#### QATrack+ File Upload: Generic Binary File

---

**Todo:** Text File upload docs

---



## DEVELOPMENT NOTES

### 6.1 Runing tests

The tests can be run by running:

```
py.test
```

in the root qcpump directory.

### 6.2 Release Checklist

- Tests all passing
- docs/release\_notes.rst updated
- Version updated in:
  - setup.py
  - qcpump/settings.py
  - qcpump-installer.iss
  - docs/install.rst (link to installer)
- **Exe built (*pyinstaller qcpump.spec*)**
  - New data files added to qcpump.spec::data\_files
  - New dependencies added to qcpump.spec::hidden\_import
- **Installer built**
  - Installer comitted to repository
- Installer tested
- Release tagged *git tag -a vX.X.X -m vX.X.X*
- Push to master with tags *git push origin master --tags*

## 6.3 Building an exe on Windows

In order to create an executable version of QCPump on Windows you need to install pyinstaller:

```
pip install pyinstaller
```

Then run:

```
pyinstaller qcpump.spec
```

## 6.4 Building the Installer

- Install Inno Setup Compiler
- Build the exe as described above
- Open qcpump-installer.iss in Inno
- Update MyAppVersion
- Build Menu -> Compile (or Ctrl+F9)
- Test that the installer works and QCPump launches after installation
- Add the installer to the repository and push.

## PUMP TYPE DEVELOPMENT


### 7.1 Tutorial

#### 7.1.1 Install QCPump

Follow the *installation instructions* and either install the precompiled executable version of QCPump using the Windows installer, or install QCPump from source

#### 7.1.2 Create a PumpType directory




Create a directory somewhere to contain your pump types. I am going to use a directory in my home directory called *Pumps*:



Now, find your QCPump settings file location file by launching QCPump and going to File->About:

Navigate to that directory and edit the settings.json file (you can just use notepad.exe or whatever your favourite editor is):

adding your new PumpType directory to the *PUMP\_DIRECTORIES* list (use forward slashes / instead of back slashes \ !):



## 7.2 Developing your own Pump Types

To create your own PumpType you will need to create a subclass of `qcpump.pumps.base.BasePump` and, at a minimum, implement a `pump` method.

First create a subdirectory in your pumps subdirectory called `HelloWorld` and in that directory create a file called `hello_world.py` with the following content:

```
from qcpump.pumps.base import BasePump


class HelloWorldPump(BasePump):

    DISPLAY_NAME = "Hello World"


    def pump(self):
        self.log_info("Hello World from your new pump!")
```

This is the simplest possible pump; all it does is log a debug message every time QCPump calls its `pump` method.

Now launch QCPump and click add a new Hello World pump:



On the `Hello World` pump page, click `Save` then `Run Pumps` (note if you have other pumps configured you may want to deactivate them first!). You should see a message logged to the status pane of your Hello World pump:



That's it! You've just written your first Pump Type. It doesn't do much, but should get you started on your journey.

## 7.3 Adding Configuration Options To Your Pump

Let's extend our Hello World Pump to log a message that can be customized by our users.

Modify your `hello_world.py` file so it has the following code:

```
from qcpump.pumps.base import BasePump, STRING

class HelloWorldPump(BasePump):

    DISPLAY_NAME = "Hello World"

    CONFIG = [
        {
            'name': 'Hello World',
            'fields': [
                {
                    'name': 'message',
                    'type': STRING,
                    'required': True,
                    'help': "Enter the message you want to see logged",
                    'default': "A default message",
                },
            ],
        },
    ]

    def pump(self):
        message = self.get_config_value("Hello World", "message")
        self.log_info(f"Hello World says: {message}")
```

By adding the `CONFIG` option, we are telling QCPump we want our users to be able to configure an option called `message` in a configuration section called `Hello World` that will have a string type (see [Config Options](#) for more option types). If you relaunch QCPump, you should now see your new configuration option. Change the default message, click *Save*, and run the pump again and you should see it log your new custom message:

## 7.4 Adding Validation To Your Pump

QCPump allows you to do some data validation to ensure things are configured correctly. Modify your `hello_world.py` to add a validation method for our message option:

```
from qcpump.pumps.base import BasePump, STRING


class HelloWorldPump(BasePump):

    DISPLAY_NAME = "Hello World"

    CONFIG = [
        {
            'name': 'Hello World',
            'validation': 'validate_message',
            'fields': [
```

(continues on next page)





(continued from previous page)

```

        {
            'name': 'message',
            'type': STRING,
            'required': True,
            'help': "Enter the message you want to see logged",
            'default': "A default message",
        },
    ],
},
]

def pump(self):
    message = self.get_config_value("Hello World", "message")
    self.log_info(f"Hello World says: {message}")

def validate_message(self, values):
    message = values['message']

    if "QCPump" not in message:
        valid = False
        message = "You must include the text 'QCPump' in your message!"
    else:
        valid = True
        message = "Thank you for including 'QCPump' in your message!"

    return valid, message

```

Here we have added a `validate_message` method that ensures our users have the text 'QCPump' in their message. The `values` variable is a dictionary with keys made up by the names of our configuration options from this section, and values made up of the current value the user has configured. Since our `Hello World` configuration section only has a single option, `values` is just a dictionary with a single key:


```

values == {
    'message': "Some configured message",
}


```

Validation methods must return an iterable of length two (e.g. a two tuple, or a list of two items) where the first item is True/False representing whether the configuration section is valid or not, and the second item is a string validation message which will be displayed to the user.

If you launch QCPump again and enter a value without the string 'QCPump' you should see your Pump will be marked as invalid:



Adding the text 'QCPump' will validate the pump:



For more information on writing QCPumps please see the information below.

## 7.5 Config Options

### 7.5.1 Option Types

You can add options of type:

**STRING** Short free form text

**BOOLEAN** A True/False drop down

**INT** An integer value

**UINT = 'uint'** A unsigned/positive integer value

**FLOAT** A decimal (floating point) value

**MULTCHOICE** A dropdown to select one option

**DIRECTORY** A path to a directory

### 7.5.2 Setting choices for Multiple Choice options

---

**Todo:** Choices Docs

---

### 7.5.3 Validation

---

**Todo:** Validation docs

---

## 7.6 Dependencies

---

**Todo:** Dependency Docs

---

## 7.7 QATrack+ Mixins

### 7.7.1 QATrackAPIMixin

---

**Todo:** QATrackAPIMixin docs

---

## 7.7.2 QATrackFetchAndPost

---

**Todo:** QATrackFetchAndPost docs

---

## 7.7.3 QATrackFetchAndPostTextFile

---

**Todo:** QATrackFetchAndPostTextFile docs

---

## 7.7.4 QATrackFetchAndPostBinarFile

---

**Todo:** QATrackFetchAndPostBinaryFile docs

---

## QC PUMP OVERVIEW




Fig. 1: QCPump Interface

QCPump is an extensible utility for extracting data from various sources (databases, filesystem etc) and moving it to other locations.

QCPump allows you to define one or more *Pumps* of various types (known as *Pump Types* e.g. *File Mover* and *Firebird DQA3*) to retrieve data from databases, filesystems, and other data sources. Typically these *Pumps* run at regular interval to look for new data, and when that data is present they extract it, potentially transform it in some way, and then upload it to QATrack+, move it somewhere else on disk, leave it alone, or take some other action.

QCPump currently includes the following Pump Types:

- *Simple File Mover & File Mover* which serve as examples for writing your pump types.

- *FirebirdDQA3* & *AtlasDQA3* which can be used for retrieving Daily QA3 data from Firebird & Atlas (SQLServer) databases.

QC Pump was built primarily as a tool for uploading data to QATrack+, but there is nothing preventing you from using it for other tasks you want to perform periodically. For example, you could *write a Pump Type of your own* that ran once a day to generate a backup file of a database. If you think your *Pump Type* would be valuable to others please contribute it back to QCPump!

## 8.1 QC Pump License

QC Pump is licensed under the [MIT License](#) and all code contributed to the QCPump project will fall under the same license.

## INDICES AND TABLES

- genindex
- modindex
- search